

anti- EDNRB antibody

Product Information

Catalog No.:	FNab02637
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Non-specific receptor for endothelin 1, 2, and 3. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

Immunogen information

Immunogen:	endothelin receptor type B
Synonyms:	ETRB
Observed MW:	53 kDa
UniprotID :	P24530

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab02637(EDNRB Antibody) at dilution of 1:50



mouse placenta tissue were subjected to SDS PAGE followed by western blot with FNab02637(EDNRB antibody) at dilution of 1:500