

## anti- EEA1 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab02639   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Protein A+G purification  |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Mouse   |
| Clonality:    | monoclonal  |
| Clone ID:     | 2H5   |
| IsoType:      | IgG1  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

Binds phospholipid vesicles containing phosphatidylinositol 3-phosphate and participates in endosomal trafficking.

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | early endosome antigen 1                                |
| Synonyms:    | early endosome antigen 1, EEA1, MST105, MSTP105, ZFYVE2 |
| Observed MW: | None  |
| Uniprot ID : | Q15075  |

### Application

|                          |            |
|--------------------------|------------|
| Reactivity:              | Human      |
| Tested Application:      | ELISA, IF  |
| Recommended dilution:IF: | 1:20-1:200 |
| Image:                   |            |

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)