

anti- ESRRB antibody

Product Information

Catalog No.:	FNab02871
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Nuclear receptor, may regulate ESR1 transcriptional activity. Induces the expression of PERM1 in the skeletal muscle.

Immunogen information

Immunogen:	estrogen-related receptor beta
Synonyms:	ERRB2, ESRL2, NR3B2
Observed MW:	45-49 kDa
UniprotID :	O95718

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-ESRRB (IP:FNab02871, 4ug;
Detection:FNab02871 1:300) with A549 cells
lysate 1040ug.



PC-13 cells were subjected to SDS PAGE
followed by western blot with FNab02871(ESRRB
antibody) at dilution of 1:1000