

anti- ESRRG antibody

Product Information

Catalog No.:	FNab02872
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Orphan receptor that acts as transcription activator in the absence of bound ligand. Binds specifically to an estrogen response element and activates reporter genes controlled by estrogen response elements(By similarity). Induces the expression of PERM1 in the skeletal muscle.

Immunogen information

Immunogen:	estrogen-related receptor gamma
Synonyms:	ERR3, ERG2, KIAA0832, NR3B3
Observed MW:	51 kDa
Uniprot ID :	P62508

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:3000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



human spleen tissue were subjected to SDS PAGE followed by western blot with FNab02872(ESRRG antibody) at dilution of 1:1000