

anti- EXOC1 antibody

Product Information

Catalog No.:	FNab02893
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Component of the exocyst complex involved in the docking of exocytic vesicles with fusion sites on the plasma membrane.

Immunogen information

Immunogen:	exocyst complex component 1
Synonyms:	SEC3, SEC3L1
Observed MW:	102 kDa
UniprotID :	Q9NV70

Application


Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, IHC, WB, IP	
Recommended dilution: WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200		
Image:		

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human gliomas using FNab02893(EXOC1 antibody) at dilution of 1:50

IP Result of anti-EXOC1 (IP:FNab02893, 3ug; Detection:FNab02893 1:1000) with mouse brain tissue lysate 3600ug.




Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





rat brain tissue were subjected to SDS PAGE followed by western blot with FNab02893(EXOC1 antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>