

anti- Factor XII antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab02939 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Protein A+G purification |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 4E7 |
| IsoType: | IgG1 |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

None

Immunogen information

| | |
|--------------|--|
| Immunogen: | coagulation factor XII(Hageman factor) |
| Synonyms: | None |
| Observed MW: | 28 kDa |
| UniprotID : | P00748 |

Application


| | |
|-----------------------|---|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB, IF, IHC |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:500 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab02939(F12 Antibody) at dilution of 1:200



K-562 cells were subjected to SDS PAGE followed by western blot with FNab02939(F12 antibody) at dilution of 1:1000