

## anti- FAM107A antibody

### Product Information

Catalog No.:	FNab02953
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

When transfected into cell lines in which it is not expressed, suppresses cell growth. May play a role in tumor development.

### Immunogen information

Immunogen:	family with sequence similarity 107, member A
Synonyms:	DRR1, TU3A
Observed MW:	None
UniprotID :	O95990

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human lung cancer using FNab02953(FAM107A antibody) at dilution of 1:50