

anti- FAM127A antibody

Product Information

Catalog No.:	FNab02962
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	family with sequence similarity 127, member A
Synonyms:	CXX1, MAR8, MAR8C, MART8
Observed MW:	40 kDa
UniprotID :	A6ZKI3

Application


Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab02962(FAM127A antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02962(FAM127A antibody) at dilution of 1:300