

anti- FAM3C antibody

Product Information

Catalog No.:	FNab02981
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in retinal laminar formation. Promotes epithelial to mesenchymal transition.

Immunogen information

Immunogen:	family with sequence similarity 3, member C
Synonyms:	FAM3C, Protein FAM3C, GS3786, Interleukin like EMT inducer, ILEI
Observed MW:	25kd
Uniprot ID :	Q92520

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas using FNab02981(FAM3C antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02981(FAM3C antibody) at dilution of 1:500