

## anti- FAM50A antibody

### Product Information

Catalog No.:	FNab02988
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May be a DNA-binding protein or transcriptional factor.

### Immunogen information

Immunogen:	family with sequence similarity 50, member A
Synonyms:	DXS9928E, HXC26, XAP5
Observed MW:	37-40 kDa
UniprotID :	Q14320

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02988(FAM50A antibody) at dilution of 1:1500