

anti- FEM1A antibody

Product Information

Catalog No.:	FNab03071
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Probable component of an E3 ubiquitin-protein ligase complex, in which it may act as a substrate recognition subunit(By similarity). May participate in antiinflammatory signaling via its interaction with PTGER4.

Immunogen information

Immunogen:	fem-1 homolog a(C. elegans)
Synonyms:	EPRAP
Observed MW:	74 kDa
UniprotID :	Q9BSK4

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab03071(FEM1A antibody) at dilution of 1:50



U-937 cells were subjected to SDS PAGE followed by western blot with FNab03071(FEM1A antibody) at dilution of 1:300