

## anti- FERMT3 antibody

### Product Information

Catalog No.:	FNab03078
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Plays a central role in cell adhesion in hematopoietic cells(PubMed:19234463, PubMed:26359933). Acts by activating the integrin beta-1-3(ITGB1, ITGB2 and ITGB3)(By similarity). Required for integrin-mediated platelet adhesion and leukocyte adhesion to endothelial cells(PubMed:19234460). Required for activation of integrin beta-2(ITGB2) in polymorphonucleargranulocytes(PMNs)(By similarity). Isoform 2 may act as a repressor of NF-kappa-B and apoptosis.

### Immunogen information

Immunogen:	fermitin family homolog 3(Drosophila)
Synonyms:	KIND3, MIG2B, URP2
Observed MW:	72-76kd
UniprotID :	Q86UX7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded human lymphoma using FNab03078(FERMT3 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03078(FERMT3 antibody) at dilution of 1:300

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)