

anti- FGFBP1 antibody

Product Information

Catalog No.:	FNab03096
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Acts as a carrier protein that release fibroblast-binding factors(FGFs) from the extracellular matrix(EM) storage and thus enhance the mitogenic activity of FGFs. Enhances FGF2 signaling during tissue repair, angiogenesis and in tumor growth.

Immunogen information

Immunogen:	fibroblast growth factor binding protein 1
Synonyms:	17 kDa HBGF binding protein, FGF binding protein 1, FGF BP, FGF BP1, FGFBP, FGFBP 1, FGFBP1, HBP17, Heparin-binding protein 17
Observed MW:	26 kDa
Uniprot ID :	Q14512

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:50-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human colon cancer using FNab03096(FGFBP1 antibody) at dilution of 1:50