

## anti- Frizzled 7 antibody

### Product Information

Catalog No.:	FNab03221
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	frizzled homolog 7(Drosophila)
Synonyms:	None
Observed MW:	64 kDa
UniprotID :	O75084

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-Frizzled 7 (IP:FNab03221, 4ug;  
Detection:FNab03221 1:600) with mouse skeletal  
muscle tissue lysate 2200ug.



mouse skeletal muscle tissue were subjected to  
SDS PAGE followed by western blot with  
FNab03221(Frizzled 7 antibody) at dilution of  
1:600