

anti- FSTL1 antibody

Product Information

Catalog No.:	FNab03232
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May modulate the action of some growth factors on cell proliferation and differentiation. Binds heparin(By similarity).

Immunogen information

Immunogen:	follistatin-like 1
Synonyms:	FRP
Observed MW:	35-45 kDa
Uniprot ID :	Q12841

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab03232(FSTL1 Antibody) at dilution of 1:50



A2780 cells were subjected to SDS PAGE followed by western blot with FNab03232(FSTL1 antibody) at dilution of 1:300