

anti- FTCD antibody

Product Information

Catalog No.: FNab03233

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Folate-dependent enzyme, that displays both transferase and deaminase activity. Serves to channel one-carbon units from formiminoglutamate to the folate pool. Binds and promotes bundling of vimentin filaments originating from the Golgi.

Immunogen information

Immunogen: formiminotransferasecyclodeaminase

Synonyms: None
Observed MW: 53kd
UniprotID: 095954

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human kidney using FNab03233(FTCD antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab03233(FTCD antibody) at dilution of 1:1000