
anti- FZR1 antibody

Product Information

Catalog No.:	FNab03265
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	fizzy/cell division cycle 20 related 1 (Drosophila)
Synonyms:	CDH1, FYR, FZR, KIAA1242
Observed MW:	55 kDa
Uniprot ID :	Q9UM11

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:20 - 1:50
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney cancer using FNa03265(FZR1 Antibody) at dilution of 1:100



K-562 cells were subjected to SDS PAGE followed by western blot with FNa03265(FZR1 antibody) at dilution of 1:1000