

anti- GABPA antibody

Product Information

Catalog No.:	FNab03280
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcription factor capable of interacting with purine rich repeats(GA repeats). Necessary for the expression of the Adenovirus E4 gene.

Immunogen information

Immunogen:	GA binding protein transcription factor, alpha subunit 60kDa
Synonyms:	E4TF1A
Observed MW:	56-60 kDa
UniprotID :	Q06546

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer using FNab03280(GABPA antibody) at dilution of 1:50



IP Result of anti-GABPA (IP:FNab03280, 4ug; Detection:FNab03280 1:1000) with HeLa cells lysate 1200ug.



K-562 cells were subjected to SDS PAGE followed by western blot with FNab03280(GABPA antibody) at dilution of 1:1000