

anti- GAD2 antibody

Product Information

Catalog No.:	FNab03292
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Catalyzes the production of GABA.

Immunogen information

Immunogen:	glutamate decarboxylase 2(pancreatic islets and brain, 65kDa)
Synonyms:	GAD65
Observed MW:	65 kDa
UniprotID :	Q05329

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-GAD65 (IP:FNab03292, 4ug;
Detection:FNab03292 1:1000) with mouse brain
tissue lysate 3440ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with FNab03292(GAD2
antibody) at dilution of 1:1000