

# anti- GATA3 antibody

#### **Product Information**

Catalog No.: FNab03361

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Background**

Transcriptional activator which binds to the enhancer of the T-cell receptor alpha and delta genes. Binds to the consensus sequence 5'-AGATAG-3'. Required for the T-helper 2(Th2) differentiation process following immune and inflammatory responses.

## Immunogen information

Immunogen: GATA binding protein 3

Synonyms: None

Observed MW: 45-50 kDa UniprotID: P23771

### **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab03361(GATA3 Antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab03361(GATA3 antibody) at dilution of 1:500