

anti- Giantin antibody

Product Information

Catalog No.:	FNab03453
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	golgin B1, golgi integral membrane protein
Synonyms:	GCP, GCP372, giantin, GOLGB1, Golgin subfamily B member 1, GOLIM1, Macrogolgin
Observed MW:	376 kDa
Uniprot ID :	Q14789

Application

Reactivity:	Human
Tested Application:	ELISA, IF
Recommended dilution:IF:	1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunofluorescent analysis of HepG2 cells using FNab03453(Giantin antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)