

---

## anti- HBZ-Specific antibody

### Product Information

Catalog No.:	FNab03780
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	hemoglobin, zeta
Synonyms:	None
Observed MW:	42kd
UniprotID :	P02008

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:4000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03780(HBZ-Specific antibody) at dilution of 1:1000