

## anti- HESX1 antibody

### Product Information

Catalog No.:	FNab03841
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Required for the normal development of the forebrain, eyes and other anterior structures such as the olfactory placodes and pituitary gland. Possible transcriptional repressor. Binds to the palindromic PIII sequence, 5'-AGCTTGAGTCTAATTGAATTAAGTGTAC-3'. HESX1 and PROP1 bind as heterodimers on this palindromic site, and, in vitro, HESX1 can antagonize PROP1 activation (By similarity).

### Immunogen information

Immunogen:	HESX homeobox 1
Synonyms:	HANF
Observed MW:	None
UniprotID :	Q9UBX0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution: IHC:	1:20-1:200


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Image:



Immunohistochemistry of paraffin-embedded human brain using FNab03841(HESX1 Antibody) at dilution of 1:50

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)