

## anti- HEXA antibody

### Product Information

Catalog No.:	FNab03843
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4F7
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Responsible for the degradation of GM2 gangliosides, and a variety of other molecules containing terminal N-acetyl hexosamines, in the brain and other tissues. The form B is active against certain oligosaccharides. The form S has no measurable activity.

### Immunogen information

Immunogen:	hexosaminidaseA(alpha polypeptide)
Synonyms:	HEXA, Hexosaminidase subunit A, TSD
Observed MW:	60 kDa
UniprotID :	P06865

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:5000 - 1:20000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03843(HEXA antibody) at dilution of 1:10000