

## anti- HMMR-Specific antibody

### Product Information

Catalog No.:	FNab03936
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	hyaluronan-mediated motility receptor(RHAMM)
Synonyms:	None
Observed MW:	70 kDa
Uniprot ID :	O75330

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:10-1:100
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab03936( HMMR-Specific Antibody) at dilution of 1:200



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03936(HMMR-Specific antibody) at dilution of 1:800