

anti- HNRNPA3 antibody

Product Information

Catalog No.:	FNab03948
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in cytoplasmic trafficking of RNA. Binds to the cis-acting response element, A2RE. May be involved in pre-mRNA splicing.

Immunogen information

Immunogen:	heterogeneous nuclear ribonucleoprotein A3
Synonyms:	2610510D13Rik, D10S102, FBRNP, hnRNP A3, HNRNPA3, HNRPA3
Observed MW:	40 kDa
UniprotID :	P51991

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunofluorescent analysis of HeLa cells using FNab03948 (HNRNPA3 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03948(HNRNPA3 Antibody) at dilution of 1:1000