

anti- HRAS-Specific antibody

Product Information

Catalog No.:	FNab04008
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	v-Ha-ras Harvey rat sarcoma viral oncogene homolog
Synonyms:	C BAS/HAS, C H RAS, C HA RAS1, CTLO, GTPaseHRas, H Ras 1, H RASIDX, Ha Ras, HAMSIV, HRAS, HRAS1, K RAS, RASH1, Transforming protein p21
Observed MW:	21 kDa
UniprotID :	P01112

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IHC, WB, FC
Recommended dilution:	WB: 1:200-1:1000; IP:1:1000-1:6000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney slide using FNab04008(HRAS-Specific Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04008(HRAS-Specific Antibody) at dilution of 1:300