

anti- HSP70 antibody

Product Information

Catalog No.:	FNab04049
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Reakground	

Background

None

Immunogen information

Immunogen:	heat shock 70kDa protein 1B
Synonyms:	None
Observed MW:	66-70 kDa
Uniprot ID :	P0DMV8

Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB,IHC		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:500			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275


www.fn-test.com





Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab04048(HSP70 Antibody) at dilution of 1:200

HeLa cells were subjected to SDS PAGE followed by western blot with FNab04049(HSP70 Antibody) at dilution of 1:1000



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com