

# anti- HSP90 antibody

## **Product Information**

Catalog No.:	FNab04051
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

HSP90, encoded by HSP90AA1, is a constitutively and ubiquitously expressed molecular chaperone that is crucial for the stability and function of many proteins. HSP90 provides chaperoning activity for client proteins; many of them are members of oncogenic pathways, indicating its implication in tumor malignancy. HSP90 mainly resides in the cytosol, while it can also be released to the extracellular space. Secreted Hsp90 is a C-terminal truncated form. It has been reported that the level of plasma Hsp90 is positively correlated with tumor malignancy in clinical cancer patients, and can be a promising diagnostic marker for tumor malignancy in clinical application.

#### **Immunogen information**

Immunogen:	heat shock protein 90kDa alpha(cytosolic), class A member 1
Synonyms:	HSP 86, HSP86, Hsp89, HSP89A, HSP90, HSP90A, HSP90AA1, HSP90N, HSPC1, HSPCA, HSPCAL1, HSPCAL4, HSPN, LAP2
Observed MW:	90 kDa
Uniprot ID :	P07900

## Application


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com



Reactivity:Human, Mouse, RatTested Application:ELISA, IF, IHC, IP, WB, FCRecommended dilution:WB: 1:1000-1:6000 IP:1:500-1:1000 IHC: 1:20-1:200 IF: 1:10-1:100Image:



 $\begin{array}{c} \text{Control 190}\\ \text{Control 1900}\\ 250\text{kd} \rightarrow \\ 150\text{kd} \rightarrow \\ 100\text{kd} \rightarrow \\ 70\text{kd} \rightarrow \\ 50\text{kd} \rightarrow \end{array}$ 

Immunohistochemistry of paraffin-embedded human liver cancer using FNab04051(HSP90 antibody) at dilution of 1:50

IP Result of anti-HSP90 (IP: FNab04051, 4ug; Detection: FNab04051 1:500) with K-562 cells lysate 4000ug.


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





HeLa cells were subjected to SDS PAGE followed by western blot with FNab04051(HSP90 Antibody) at dilution of 1:3000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com