

anti- IA-2/PTPRN antibody

Product Information

Catalog No.:	FNab04093
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
IsoType:	IgG1
Clone ID:	8C8
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	protein tyrosine phosphatase, receptor type, N
Synonyms:	IA 2/PTP, IA2, IA-2, ICA 512, ICA3, ICA512, Islet cell antigen 512, Islet cell autoantigen 3, PTP IA 2, PTPRN, R PTP N
Observed MW:	60-65kd
Uniprot ID :	Q16849

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:200-1:800; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab04093(PTPRN Antibody) at dilution of 1:400 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



human brain tissue were subjected to SDS PAGE followed by western blot with FNab04093(PTPRN antibody) at dilution of 1:500