

## anti- IFT43 antibody

### Product Information

Catalog No.:	FNab04158
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Component of IFT complex A(IFT-A) involved in retrograde ciliary transport along microtubules from the ciliary tip to the base.

### Immunogen information

Immunogen:	chromosome 14 open reading frame 179
Synonyms:	C14orf179
Observed MW:	None
UniprotID :	Q96FT9

### Application


Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human testis slide using FNab04158( IFT43 antibody) at dilution of 1:50