

anti- IFT80 antibody

Product Information

Catalog No.:	FNab04161
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Component of the intraflagellartransport(IFT) complex B, which is essential for the development and maintenance of motile and sensory cilia.

Immunogen information

Immunogen:	intraflagellar transport 80 homolog(Chlamydomonas)
Synonyms:	KIAA1374, WDR56
Observed MW:	88-90 kDa
UniprotID :	Q9P2H3

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab04161(IFT80 Antibody) at dilution of 1:1000