

## anti- IL11 antibody

### Product Information

Catalog No.:	FNab04216
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Cytokine that stimulates the proliferation of hematopoietic stem cells and megakaryocyte progenitor cells and induces megakaryocyte maturation resulting in increased platelet production(PubMed:2145578). Also promotes the proliferation of hepatocytes in response to liver damage. Binding to its receptor formed by IL6ST and either IL11RA1 or IL11RA2 activates a signaling cascade that promotes cell proliferation(PubMed:12919066). Signaling leads to the activation of intracellular protein kinases and the phosphorylation of STAT3.

### Immunogen information

Immunogen:	interleukin 11
Synonyms:	Adipogenesis inhibitory factor, AGIF, IL 11, IL11, interleukin 11, Oprelvekin
Observed MW:	25 kDa
Uniprot ID :	P20809

### Application

Reactivity:	Human, Mouse
-------------	--------------

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:1000

Image:



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab04216(IL11 Antibody) at dilution of 1:500

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)