

## anti- INPP5A antibody

### Product Information

Catalog No.:	FNab04326
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Major isoenzyme hydrolyzing the calcium-mobilizing second messenger Ins(1,4,5)P<sub>3</sub>, this is a signal-terminating reaction.

### Immunogen information

Immunogen:	inositol polyphosphate-5-phosphatase, 40kDa
Synonyms:	None
Observed MW:	43-48kd
UniprotID :	Q14642

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04326(INPP5A antibody) at dilution of 1:500