

## anti- IRF6 antibody

### Product Information

Catalog No.:	FNab04393
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Probable DNA-binding transcriptional activator. Key determinant of the keratinocyte proliferation-differentiation switch involved in appropriate epidermal development(By similarity). Plays a role in regulating mammary epithelial cell proliferation(By similarity). May regulate WDR65 transcription(By similarity).

### Immunogen information

Immunogen:	interferon regulatory factor 6
Synonyms:	IRF6,LPS,OFC6,PIT,PPS,PPS1,VWS,VWS1
Observed MW:	55-65 kDa
Uniprot ID :	O14896

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)

Image:



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab04393( IRF6 Antibody) at dilution of 1:600