

## anti- KIF7 antibody

### Product Information

Catalog No.:	FNab04570
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Essential for hedgehog signaling regulation: acts as both a negative and positive regulator of sonic hedgehog(Shh) and Indian hedgehog(Ihh) pathways, acting downstream of SMO, through both SUFU-dependent and-independent mechanisms(PubMed:21633164). Involved in the regulation of microtubular dynamics. Required for proper organization of the ciliary tip and control of ciliary localization of SUFU-GLI2 complexes(By similarity). Required for localization of GLI3 to cilia in response to Shh. Negatively regulates Shh signaling by preventing inappropriate activation of the transcriptional activator GLI2 in the absence of ligand. Positively regulates Shh signaling by preventing the processing of the transcription factor GLI3 into its repressor form. In keratinocytes, promotes the dissociation of SUFU-GLI2 complexes, GLI2 nuclear translocation and Shh signaling activation(By similarity). Involved in the regulation of epidermal differentiation and chondrocyte development(By similarity).

### Immunogen information

Immunogen:	kinesin family member 7
Synonyms:	KIF7, kinesin family member 7, Kinesin like protein KIF7, UNQ340
Observed MW:	150 kDa
UniprotID :	Q2M1P5

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)


## Application

Reactivity: Human


Tested Application: ELISA, IP, WB

Recommended dilution: WB: 1:500-1:2000; IP: 1:500-1:1000

Image:



IP Result of anti-KIF7 (IP: FNab04570, 4ug;  
Detection: FNab04570 1:500) with HEK-293 cells  
lysate 2000ug.



HEK-293 cells were subjected to SDS PAGE  
followed by western blot with FNab04570( KIF7  
Antibody) at dilution of 1:1000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)