

---

## anti- KRT23 antibody

### Product Information

Catalog No.:	FNab04651
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	keratin 23(histone deacetylase inducible)
Synonyms:	CK 23, CK23, Cytokeratin 23, HAIK1, K23, Keratin 23, KRT23
Observed MW:	48 kDa
Uniprot ID :	Q9C075

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab04651(KRT23 Antibody) at dilution of 1:50



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab04651(KRT23 antibody) at dilution of 1:1000