

anti- KRT40 antibody

Product Information

Catalog No.:	FNab04655
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in late hair differentiation.

Immunogen information

Immunogen:	keratin 40
Synonyms:	CK 40, Cytokeratin 40, FLJ36600, K40, KA36, keratin 40, KRT40, Type I hair keratin Ka36
Observed MW:	48 kDa
Uniprot ID :	Q6A162

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



A431 cells were subjected to SDS PAGE followed by western blot with FNab04655(KRT40 Antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com