

## anti- LAYN antibody

### Product Information

Catalog No.:	FNab04711
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Receptor for hyaluronate.

### Immunogen information

Immunogen:	layilin
Synonyms:	layilin, LAYN
Observed MW:	43kd
UniprotID :	Q6UX15

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:10-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunofluorescent analysis of HepG2 cells, using FNab04711 (LAYN antibody) at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red)



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04711(LAYN antibody) at dilution of 1:1000