

## anti- LILRA2 antibody

### Product Information

Catalog No.:	FNab04775
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Does not bind class I MHC antigens.

### Immunogen information

Immunogen:	leukocyte immunoglobulin-like receptor, subfamily A(with TM domain), member 2
Synonyms:	ILT1, LIR7
Observed MW:	53kd
UniprotID :	Q8N149

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab04775(LILRA2 Antibody) at dilution of 1:200



human liver tissue were subjected to SDS PAGE followed by western blot with FNab04775(LILRA2antibody) at dilution of 1:500