

# anti- liver Arginase antibody

#### **Product Information**

Catalog No.: FNab04800

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Arginase-1(Liver arginase) belongs to the arginase family. ARG1 is a novel immunohistochemical marker of hepatocellular differentiation in fine needle aspiration cytology and a marker of hepatocytes and hepatocellular neoplasms. ARG1 is closely associated with alternative macrophage activation(PMID:12098359) and ARG1 has been shown to protectmotor neurons from trophic factor deprivation and allow sensory neurons to overcome neurite outgrowth inhibition by myelin proteins(PMID: 20071539). It can exsit as a homotrimer(PMID:16141327) and it has 3 isoforms produced by alternative splicing. Defects in ARG1 are the cause of argininemia(ARGIN). Deletion or TNF-mediated restriction of ARG1 unleashes the production of NO by NOS2, which is critical for pathogen control. (PMID:27117406).

## **Immunogen information**

Immunogen: arginase, liver

Synonyms: ARG1, Arginase 1, liver Arginase, Liver type arginase

Observed MW: 35 kDa UniprotID: P05089

### 1

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>



# **Application**


Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC, IP, IF, FC, WB


Recommended dilution: WB: 1:1000-1:10000; IP: 1:500-1:1000; IHC: 1:100-1:400; IF: 1:20-

1:200

Image:




Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab04800(ARG1 Antibody) at dilution of 1:200




Immunofluorescent analysis of HepG2 cells using FNab04800( ARG1 Antibody) at dilution of 1:50 and Alexa Fluor 488-congugated Goat Anti-Rabbit IgG(H+L)





IP Result of anti-ARG1 (IP:FNab04800, 4ug; Detection:FNab04800 1:700) with mouse liver tissue lysate 5200ug.



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab04800(ARG1 antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com