

anti- LSM5 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab04876 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Plays a role in U6 snRNP assembly and function. Binds to the 3' end of U6 snRNA, thereby facilitating formation of the spliceosomal U4/U6 duplex formation in vitro.

Immunogen information

| | |
|--------------|--|
| Immunogen: | LSM5 homolog, U6 small nuclear RNA associated(S. cerevisiae) |
| Synonyms: | None |
| Observed MW: | 10 kDa |
| UniprotID : | Q9Y4Y9 |

Application


| | |
|---------------------------|--------------|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04876(LSM5 antibody) at dilution of 1:600