

## anti- LSR antibody

### Product Information

Catalog No.:	FNab04880
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Probable role in the clearance of triglyceride-rich lipoprotein from blood. Binds chylomicrons, LDL and VLDL in presence of free fatty acids and allows their subsequent uptake in the cells(By similarity).

### Immunogen information

Immunogen:	lipolysis stimulated lipoprotein receptor
Synonyms:	LISCH
Observed MW:	56 kDa, 68 kDa
UniprotID :	Q86X29

### Application


Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Rat liver tissue were subjected to SDS PAGE followed by western blot with FNab04880(LSR antibody) at dilution of 1:300