

anti- LTBP1 antibody

Product Information

Catalog No.:	FNab04885
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in the assembly, secretion and targeting of TGFβ1 to sites at which it is stored and/or activated. May play critical roles in controlling and directing the activity of TGFβ1. May have a structural role in the extra cellular matrix(ECM).

Immunogen information

Immunogen:	latent transforming growth factor beta binding protein 1
Synonyms:	LTBP 1, LTBP1, TGF beta1 BP 1
Observed MW:	150 kDa
UniprotID :	Q14766

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Fnab04885(LTBP1 antibody) at dilution of 1:200



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab04885(LTBP1 Antibody) at dilution of 1:2000