

## anti- MAGED4 antibody

### Product Information

Catalog No.:	FNab04949
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	melanoma antigen family D, 4
Synonyms:	KIAA1859, MAGE D4 antigen, MAGE E1, MAGE E1 antigen, MAGE1, MAGED4, MAGED4A, MAGED4B, MAGEE1, melanoma antigen family D, 4, Melanoma associated antigen D4
Observed MW:	42-60 kDa
Uniprot ID :	Q96JG8

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:5000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab04949(MAGED4 antibody) at dilution of 1:5000