

anti- MAK16 antibody

Product Information

Catalog No.:	FNab04957
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	MAK16 homolog(S. cerevisiae)
Synonyms:	MAK16, MAK16 homolog (S. cerevisiae), MAK16L, NNP78, Protein MAK16 homolog, Protein RBM13, RBM13
Observed MW:	48 kDa
UniprotID :	Q9BXY0

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:4000; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using FNab04957 (MAK16 antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04957(MAK16 antibody) at dilution of 1:500