

anti- MARCH5 antibody

Product Information

Catalog No.:	FNab05003
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Mitochondrial E3 ubiquitin-protein ligase that plays a crucial role in the control of mitochondrial morphology by acting as a positive regulator of mitochondrial fission. May play a role in the prevention of cell senescence acting as a regulator of mitochondrial quality control. Promotes ubiquitination of FIS1, DNM1L and MFN1.

Immunogen information

Immunogen:	membrane-associated ring finger(C3HC4) 5
Synonyms:	RNF153
Observed MW:	31 kDa/65-70 kDa
Uniprot ID :	Q9NX47

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using FNab05003(MARCH5 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05003(MARCH5 antibody) at dilution of 1:500