

# anti- MCSF antibody

#### **Product Information**

Catalog No.:	FNab05068
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

# Background

None

## **Immunogen information**

Immunogen:	colony stimulating factor 1(macrophage)
Synonyms:	None
Observed MW:	70 kDa
UniprotID :	P09603

### Application


Reactivity:	Human		
Tested Application:	ELISA, WB		
Recommended dilution:WB: 1:500-1:2000			
Image:			

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





U-937 cells were subjected to SDS PAGE followed by western blot with FNab05068(MCSF antibody) at dilution of 1:200

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>