

## anti- MFAP5 antibody

### Product Information

Catalog No.:	FNab05147
Size:	100µ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a role in hematopoiesis. In the cardiovascular system, could regulate growth factors or participate in cell signaling in maintaining large vessel integrity(By similarity).

### Immunogen information

Immunogen:	microfibrillar associated protein 5
Synonyms:	MAGP 2, MAGP2, MFAP 5, MFAP5, MP25
Calculated MW:	20kDa
Uniprot ID :	Q13361

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:1000-1:3000; IHC: 1:200-1:1000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human placenta using FNab05147(MFAP5 antibody) at dilution of 1:100



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab05147(MFAP5 antibody) at dilution of 1:600