

## anti- MGP antibody

### Product Information

Catalog No.:	FNab05171
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8C3
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Associates with the organic matrix of bone and cartilage. Thought to act as an inhibitor of bone formation.

### Immunogen information

Immunogen:	matrix Gla protein
Synonyms:	GIG36, matrix Gla protein, MGLAP, MGP, NTI
Observed MW:	13 kDa
Uniprot ID :	P08493

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:500-1:1500
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human lung using FNab05171(MGP antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab05171(MGP antibody) at dilution of 1:1000